

# RESPIRATORY DEVICES

Product	Valve	Mask*	Flow Signal	Chamber Volume	Notes	Company
<b>HAND-HELD SPACERS</b>						
<b>Microchamber</b>	none	none	none	60mL	<ul style="list-style-type: none"> <li>portable and easy to carry</li> <li>latex-free</li> </ul>	<i>Respiratory Delivery Systems</i> (978) 970-1947 rdsusa.com
<b>VALVED HOLDING CHAMBER</b>						
<b>AeroChamber Plus Flow-Vu</b>	one-way	none	audio + visual	149mL	<ul style="list-style-type: none"> <li>delivers more medication than original AeroChamber</li> <li>anti-static valved holding chamber</li> <li>low effort valves</li> <li>latex-free</li> </ul>	<i>Allergan</i> (800) 678-1605 aerochambervhc.com
		small	visual			
		medium				
<b>InspiraChamber</b>	one-way	none	audio	—	<ul style="list-style-type: none"> <li>available with SootherMask (pacifier-ready) and InspiraMask</li> <li>sensitive audio and visual indicator</li> <li>latex-free</li> <li>dishwasher safe</li> </ul>	<i>Lupin</i> (800) 399-2561 inspirachamber.com
		small	audio + visual			
		medium				
<b>OptiChamber Diamond</b>	one-way	none	audio	140mL	<ul style="list-style-type: none"> <li>visible expiratory valve</li> <li>detachable LiteTouch VHC mask</li> <li>flat bottom to prevent rolling</li> </ul>	<i>Philips Respironics</i> (800) 345-6443 usa.philips.com
		small				
		medium				
<b>Vortex</b>	one-way	none	visual	194mL	<ul style="list-style-type: none"> <li>non-electrostatic metal chamber improves dose delivery consistency</li> <li>latex-free</li> <li>dishwasher safe</li> <li>compact</li> <li>universal receptacle</li> </ul>	<i>Pari</i> (800) 327-8632 pari.com
		infant				
		pediatric				
		adult, regular				
<b>Flexichamber</b>	one-way	none	visual	365cm <sup>3</sup>	<ul style="list-style-type: none"> <li>do not use abrasive cleansers or dishwasher</li> <li>bag is reusable</li> <li>collapsible and compact</li> <li>latex-free</li> </ul>	<i>Cerecor</i> (410) 522-8707 cerecor.com
		child, small				
		child, large				
		adult, small				

## NOTES

\*none = available without mask.

Suggested ages given; all masks should be measured according to manufacturer's recommendations and fit to form a tight seal with face. Not an inclusive list of products or medical devices. Please contact company for full product information.